

Translated English of Chinese Standard: GB/T35829-2018

www.ChineseStandard.net → Buy True-PDF → Auto-delivery.

Sales@ChineseStandard.net

GB

NATIONAL STANDARD OF THE
PEOPLE'S REPUBLIC OF CHINA

ICS 71.100.70

Y 42

GB/T 35829-2018

**Determination of Four Naphthalenediols in Cosmetics
- High Performance Liquid Chromatography**

化妆品中 4 种萘二酚的测定 高效液相色谱法

Issued on: February 06, 2018

Implemented on: September 01, 2018

**Issued by: General Administration of Quality Supervision, Inspection and
Quarantine;**

Standardization Administration of PRC.

Table of Contents

Foreword.....3

Introduction4

1 Scope.....5

2 Normative References5

3 Principles5

4 Reagents and Materials6

5 Apparatus.....6

6 Analytical Procedures7

7 Results Calculation8

8 Recovery Rate and Precision9

9 Allowable Deviation.....9

Appendix A (Informative) Information on Naphthalenediol Standard Substance
..... 10

Appendix B (Informative) Liquid Chromatogram and Spectrogram of Four
Naphthalenediols Standard Substance Solutions 11

Appendix C (Informative) Confirmation Test 14

Determination of Four Naphthalenediols in Cosmetics

- High Performance Liquid Chromatography

1 Scope

This Standard specifies high performance liquid chromatography for the determination of 1,7-naphthalenediol, 2,3-naphthalenediol, 1,5-naphthalenediol and 2,7-naphthalenediol in cosmetics.

This Standard is applicable to the determination of 1,7-naphthalenediol, 2,3-naphthalenediol, 1,5-naphthalenediol and 2,7-naphthalenediol in non-wax-based cream, lotion, gel, water and powder cosmetics.

The method detection limit of this Standard: the method detection limit of four naphthalenediols is 2.5mg/kg.

The method quantification limit of this Standard: the method quantification limit of four naphthalenediols is 7.5mg/kg.

2 Normative References

The following documents are essential to the application of this document. For the dated documents, only the versions with the dates indicated are applicable to this document; for the undated documents, only the latest version (including all the amendments) are applicable to this document.

GB/T 6682 Water for Analytical Laboratory Use - Specification and Test Methods

3 Principles

The samples of cream, lotion, gel and water cosmetics shall be directly extracted by 95% ethanol; the samples of powder cosmetics shall be extracted by the mixture of 95% ethanol and 0.1% acetic acid solutions; the extracting solution shall be determined by high-performance liquid chromatography after the filtering and centrifuging. According to the retention time and spectrogram qualitative, use the external standard quantitative method; if necessary, use liquid chromatography-mass spectrometry/mass spectrometry to verify.

6 Analytical Procedures

6.1 Specimen processing

6.1.1 Cream, lotion and water specimens

Take 0.5g (accurate to 0.001g) of specimen; place it into 10mL stoppered colorimetric tube; add 4mL of 95% ethanol (4.3); shake and mix on the vortex oscillator (5.3); then perform ultrasonic extraction for 15min; cool off to the room temperature; use 95% ethanol (4.3) to make constant volume to 5mL. Take 2mL of solution into the centrifuge tube; centrifuge for 2min at 10000r/min; filter the supernatant through a filter membrane (4.8) then determined by the high-performance liquid chromatography (5.1).

6.1.2 Powder specimen

Take 0.5g (accurate to 0.001g) of specimen; place it into 10mL stoppered colorimetric tube; add 2mL of 95% ethanol (4.3); shake and mix on the vortex oscillator (5.3); then add 2mL of 0.1% acetic acid solution (4.4); perform ultrasonic extraction for 15min; cool off to the room temperature; use 95% ethanol (4.3) to make constant volume to 5mL. Take 2mL of solution into the centrifuge tube; centrifuge for 2min at 10000r/min; filter the supernatant through a filter membrane (4.8) then determined by the high-performance liquid chromatography (5.1).

6.2 Determination

6.2.1 Reference operating conditions of high-performance liquid chromatograph

The reference operating conditions of high-performance liquid chromatograph are as follows:

- a) Chromatographic column: C₁₈ column, 5μm, 250mm×4.6mm (inner diameter), or the equivalent;
- b) Mobile phase: methanol (4.1) and 0.1% acetic acid solution (4.4), high-performance liquid chromatograph gradient elution conditions are shown in Table 1;
- c) Flow rate: 1.0mL/min;
- d) Column temperature: 30°C;
- e) Wavelength: 230nm;
- f) Sampling-injecting volume: 20μL.

$$X_i = \frac{(c_i - c_0) \times V \times 1\,000}{m \times 1\,000} \times f$$

Where:

X_i – the content of the test composition in the specimen; in mg/kg;

c_i – the mass concentration of the test composition in the specimen solution obtained from the standard curve; in $\mu\text{g/mL}$;

c_0 – the mass concentration of the test composition in the blank sample obtained from the standard curve; in $\mu\text{g/mL}$;

V – the constant volume of the specimen solution, in mL;

m – specimen mass, in g;

f – dilution factor.

The results shall be retained two significant figures.

8 Recovery Rate and Precision

In the range of add concentration of 7.5mg/kg ~75mg/kg, the recovery rate is 80%~110%; the relative standard deviation is less than 10%.

9 Allowable Deviation

The absolute difference between two independent determinations obtained under repeatability conditions shall not exceed 10% of the arithmetic mean.

Appendix C

(Informative)

Confirmation Test

C.1 Liquid chromatography reference conditions

The liquid chromatography reference conditions are as follows:

- a) Chromatographic column: Luna C₁₈(2), 150×2.0mm (inner diameter), 3μm; or the equivalent;
- b) Mobile phase: A – methanol; B – water; gradient elution: 0.0min~10.0min, 25%~65%A; 10.1min~30.0min, 65%A.
- c) Flow rate: 150μL/min;
- d) Column temperature: 25°C;
- e) Sample-injecting volume: 10μL.

C.2 Mass spectrometry reference conditions

The mass spectrometry reference conditions are as follows:

- a) Ion source: electrospray ion source;
- b) Ionization mode: negative ion mode;
- c) Ionization voltage: -4500V;
- d) Atomizing gas pressure: 0.345MPa (50psi);
- e) Air curtain gas pressure: 0.241MPa (35psi);
- f) Auxiliary heating gas pressure: 0.483MPa (70psi);
- g) Auxiliary gas temperature: 500°C;
- h) Detection mode: Multiple Reaction Monitoring (MRM) mode. The mass spectrometric analysis parameters of four naphthalenediols can refer to Table C.1; the multiple reaction monitoring mass chromatograms of four naphthalenediols can refer to Figure C.1.

C.3 Qualitative determination

This is an excerpt of the PDF (Some pages are marked off intentionally)

Full-copy PDF can be purchased from 1 of 2 websites:

1. <https://www.ChineseStandard.us>

- SEARCH the standard ID, such as GB 4943.1-2022.
- Select your country (currency), for example: USA (USD); Germany (Euro).
- Full-copy of PDF (text-editable, true-PDF) can be downloaded in 9 seconds.
- Tax invoice can be downloaded in 9 seconds.
- Receiving emails in 9 seconds (with download links).

2. <https://www.ChineseStandard.net>

- SEARCH the standard ID, such as GB 4943.1-2022.
- Add to cart. Only accept USD (other currencies - <https://www.ChineseStandard.us>).
- Full-copy of PDF (text-editable, true-PDF) can be downloaded in 9 seconds.
- Receiving emails in 9 seconds (with PDFs attached, invoice and download links).

Translated by: Field Test Asia Pte. Ltd. (Incorporated & taxed in Singapore. Tax ID: 201302277C)

About Us (Goodwill, Policies, Fair Trading...): <https://www.chinesestandard.net/AboutUs.aspx>

Contact: Wayne Zheng, Sales@ChineseStandard.net

Linkin: <https://www.linkedin.com/in/waynezhengwenrui/>

----- The End -----