Translated English of Chinese Standard: GB/T5493-2008

<u>www.ChineseStandard.net</u> → Buy True-PDF → Auto-delivery.

<u>Sales@ChineseStandard.net</u>

GB

NATIONAL STANDARD OF THE PEOPLE'S REPUBLIC OF CHINA

ICS 67.040

X 10

GB/T 5493-2008

Replacing GB/T 5493-1985

Inspection of grain and oils - Determination of type purity and their mixture

粮油检验 类型及互混检验

Issued on: November 04, 2008 Implemented on: January 20, 2009

Issued by: General Administration of Quality Supervision, Inspection and Quarantine;
Standardization Administration of the People's Republic of China.

Table of Contents

Foreword	. 3
1 Scope	.4
2 Normative references	.4
3 Reagents	. 4
4 Instruments and utensils	.4
5 Lighting requirements	.5
6 Operation procedures	.5
7 Result calculations	.6

Inspection of grain and oils - Determination of type purity and their mixture

1 Scope

This Standard specifies the reagents, instruments and utensils, lighting requirements, operating procedures and result calculations for grain and oil varieties and their mixture.

This Standard applies to the determination of mixture of grains and oils of the same type and different grain forms, grain natures, and grain colors, as well as the inspection of mixture of different types of grain and oils.

2 Normative references

The terms in the following documents become the terms of this Standard by reference to this Standard. For dated references, all subsequent amendments (not including errata content) or revisions do not apply to this standard. However, parties to agreements that are based on this Standard are encouraged to study whether the latest versions of these documents can be used. For undated references, the latest edition applies to this Standard.

GB 1350, Paddy

GB 1351, Wheat

GB 1354, Milled rice

GB/T 5494, Inspection of grain and oils - Determination of foreign matter and unsound kernels of grain and oilseeds

GB/T 22505, Inspection of grain and oils - Sensory analysis environment - Lighting

3 Reagents

0.1% iodine-ethanol solution.

4 Instruments and utensils

- **4.1** Electronic balance: division value 0.01 g.
- **4.2** Experimental hulling machine.

6.3 Dyeing inspection of glutinous rice and non-glutinous rice

After the paddy is husked and whitened, make it into rice samples conforming to the national standard level 3 (GB 1354); take 200 whole grains out indiscriminately; wash with clean water; soak in 0.1% iodine-ethanol solution for about 1 min; then, wash with distilled water and observe the coloring of the rice grains. Glutinous rice grains are brownish red; non-glutinous rice grains are blue. Pick out the number of grains mixed with different types (n).

7 Result calculations

7.1 Calculation of mixing ratio of indica rice, japonica rice and glutinous rice

Calculate the mixing ratio of indica rice, japonica rice and glutinous rice according to Formula (1):

$$X = \frac{m_1}{m_a} \times 100 \qquad \qquad \dots \tag{1}$$

Where:

X – mixing ratio, that is, the mass fraction of mixed paddy of different types in the total sample, %;

 m_1 – the mass of paddy of different types, in grams (g);

ma – sample mass, in grams (g).

Under repeatability conditions, the absolute difference of the two independent test results obtained is not more than 1%, and the average is the test result. Retain test results to integer digits.

7.2 Calculation of the ratio of different color kernel

Calculate the ratio of different color kernel according to Formula (2):

$$Y = \frac{m_2}{m_b} \times 100 \qquad \qquad \cdots \qquad (2)$$

Where:

Y – ratio of different color kernel, that is, the mass fraction of different color kernel in the total sample, %;

m₂ – mass of different color kernel, in grams (g);

 m_b – sample mass, in grams (g).

This is an excerpt of the PDF (Some pages are marked off intentionally)

Full-copy PDF can be purchased from 1 of 2 websites:

1. https://www.ChineseStandard.us

- SEARCH the standard ID, such as GB 4943.1-2022.
- Select your country (currency), for example: USA (USD); Germany (Euro).
- Full-copy of PDF (text-editable, true-PDF) can be downloaded in 9 seconds.
- Tax invoice can be downloaded in 9 seconds.
- Receiving emails in 9 seconds (with download links).

2. https://www.ChineseStandard.net

- SEARCH the standard ID, such as GB 4943.1-2022.
- Add to cart. Only accept USD (other currencies https://www.ChineseStandard.us).
- Full-copy of PDF (text-editable, true-PDF) can be downloaded in 9 seconds.
- Receiving emails in 9 seconds (with PDFs attached, invoice and download links).

Translated by: Field Test Asia Pte. Ltd. (Incorporated & taxed in Singapore. Tax ID: 201302277C)

About Us (Goodwill, Policies, Fair Trading...): https://www.chinesestandard.net/AboutUs.aspx

Contact: Wayne Zheng, Sales@ChineseStandard.net

Linkin: https://www.linkedin.com/in/waynezhengwenrui/

---- The End -----